



OpenDiscovery

www.opendiscovery.org.uk

Version: 1.0

Contacts: gareth.price@warwick.ac.uk; a.marsh@warwick.ac.uk

Example 1: Statin + HMG-CoA Reductase

You can follow along with this example by downloading the initial PDB from [here](#). Before starting, watch the video tutorial provided by the authors of AutoDock Vina, [here](#).

Preparing the files

We first need to split the crystal structure of the complex into separate receptor and ligand files. Open up the file 1T02.pdb in your favourite text editor and look at the file. The first section are REMARKS about the structure:

```
HEADER      OXIDOREDUCTASE              07-APR-04   1T02
TITLE       CRYSTAL STRUCTURE OF A STATIN BOUND TO CLASS II HMG-COA
TITLE       2 REDUCTASE
COMPND      MOL_ID: 1;
COMPND      2 MOLECULE: 3-HYDROXY-3-METHYLGLUTARYL-COENZYME A REDUCTASE;
COMPND      3 CHAIN: A, B;
COMPND      4 SYNONYM: HMG-COA REDUCTASE;
COMPND      5 EC: 1.1.1.88;
COMPND      6 ENGINEERED: YES
SOURCE      MOL_ID: 1;
SOURCE      2 ORGANISM_SCIENTIFIC: PSEUDOMONAS MEVALONII;
SOURCE      3 ORGANISM_TAXID: 32044;
SOURCE      4 GENE: MVAA;
SOURCE      5 EXPRESSION_SYSTEM: ESCHERICHIA COLI BL21(DE3);
SOURCE      6 EXPRESSION_SYSTEM_TAXID: 469008;
SOURCE      7 EXPRESSION_SYSTEM_STRAIN: BL21(DE3);
SOURCE      8 EXPRESSION_SYSTEM_VECTOR_TYPE: PLASMID;
SOURCE      9 EXPRESSION_SYSTEM_PLASMID: PKK177-3
KEYWDS      STATIN, HMG-COA REDUCTASE, COMPLEX, OXIDOREDUCTASE
EXPDTA      X-RAY DIFFRACTION
AUTHOR      L.TABERNERO,V.W.RODWELL,C.STAUFFACHER
REVDAT      2 24-FEB-09 1T02 1 VERSN
REVDAT      1 03-AUG-04 1T02 0
JRNL        AUTH L.TABERNERO,V.W.RODWELL,C.V.STAUFFACHER
JRNL        TITL CRYSTAL STRUCTURE OF A STATIN BOUND TO A CLASS II
JRNL        TITL 2 HYDROXYMETHYLGLUTARYL-COA REDUCTASE.
JRNL        REF J.BIOL.CHEM. V. 278 19933 2003
JRNL        REFN ISSN 0021-9258
JRNL        PMID 12621048
JRNL        DOI 10.1074/JBC.M213006200
REMARK      1
REMARK      2
REMARK      2 RESOLUTION. 2.60 ANGSTROMS.
```

Delete all the header lines, up to line 765 which is the first atom of the structure.

Copy all lines with ATOM at the start into its own file, named [hmg_coa.pdb](#). Include the TER statement in line 5529. Then, copy the subsequent lines starting with HETATM and LVA in the the fourth column into its own file, [LVA.pdb](#):

Acknowledge the warnings

Save the PDBQT

Define the Grid Box and remember the size and center grid box coordinates

Open a new file and enter the box coordinates and dimensions such as:

```
center_x = 83.033
center_y = 125.978
center_z = 106.748

size_x = 16
size_y = 14
size_z = 14
```

Save it as `hmg_coa.conf`.

You should have a folder that looks like:

```

├─ hmg_coa.conf
├─ pdb
├─ └─ LVA.pdb
├─ receptor
├─ └─ hmg_coa.pdb
├─ └─ hmg_coa.pdbqt

```

Running ODScreen

We can now run the screening. First, let's check that everything is installed properly. In terminal, `cd` to the Protocol Folder (i.e. the folder with `odscreen.py`, `odparam.py` etc.) and run `python odcheck.py`. Check that there are no failures, and follow the installation guides if there are. You can ignore the warnings. Now, we can run the screening protocol. Navigate to the folder where `hmg_coa.conf` is located (from before). Now we can run `odscreen.py`. Make sure you use the correct path to the `odscreen.py` file.

```
python /protocolfolder/odscreen.py -d . -r hmg_coa -i pdb -c hmg_coa.conf
```

This is the result:

```

# ----- #
#          OPEN DISCOVERY          #
#          Screening Module         #
# ----- #
# Version: 1.0                     #
# URL:    www.opendiscovery.org.uk #
# Contacts: gareth.price@warwick.ac.uk #
#          a.marsh@warwick.ac.uk    #
# ----- #
# LigDir: /Desktop/example-statin  #
#          Receptor Name: hmg_coa   #
#          Input Type: pdb          #
#          Conf: hmg_coa.conf       #
#          Exhaustivness: 20        #
# ----- #
# Time Started: Mon, 22 Jul 2013 00:15:20 #
# ----- #

          PDB -> SMILES
Writing smiles/LVA.txt

          IMAGES
Writing images/LVA.svg

          MOL
Writing mol/LVA.mol

          MOL2
Writing mol2/LVA.mol2

          MINIMISATION
Minimising LVA

          PDBQT PREPARATION
Writing LVA.pdbqt

          SCREENING
Processing LVA

          EXTRACTING
Processing results/LVA/

          PDB -> MOL2
Writing results-mol2/LVA.mol2

```

```
SUMMARISING
Summarising LVA
```

```
MAKING COMPLEXES
Writing LVA
```

```
# ----- #
#           FINISHED           #
#       Time Taken: 231.40 seconds       #
# ----- #
```

The resulting files are:

```
.
├── combined.mol
├── hmg_coa.conf
├── images
│   └── LVA.svg
├── mol
│   └── LVA.mol
├── mol2
│   └── LVA.mol2
├── od_log.txt
├── pdb
│   └── LVA.pdb
├── pdb-minimised
│   └── LVA.pdb
├── pdbqt
│   └── LVA.pdbqt
├── receptor
│   ├── hmg_coa.pdb
│   └── hmg_coa.pdbqt
├── results
│   ├── LVA
│   │   ├── LVA_mode_1.pdb
│   │   ├── LVA_mode_10.pdb
│   │   ├── LVA_mode_2.pdb
│   │   ├── LVA_mode_3.pdb
│   │   ├── LVA_mode_4.pdb
│   │   ├── LVA_mode_5.pdb
│   │   ├── LVA_mode_6.pdb
│   │   ├── LVA_mode_7.pdb
│   │   ├── LVA_mode_8.pdb
│   │   ├── LVA_mode_9.pdb
│   │   ├── hmg_coa.pdb
│   │   ├── log.txt
│   │   └── out.pdbqt
│   ├── complexes
│   │   └── LVA.pdb
│   ├── summary-sorted.csv
│   ├── summary.csv
│   └── summary.txt
├── results-mol2
│   └── LVA.mol2
├── smiles
│   ├── LVA.txt
└── smiles.txt
```

From this point on you can analyse the results. A starting point is to look at the binding energies, sorted and summarised in the [summary-sorted.csv](#) file. This should open in any text editor or programs like Excel.